

Mitotan

Catalog No: #S1199


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mitotan
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR
Calculated MW	320.04
Formulation	C14H10Cl4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: ClC(Cl)C(C1=CC=C(Cl)C=C1)C1=CC=CC=C1Cl Pathways: Neuroscience Receptor: Adrenodoxin Binding
 pt: 398.9°C Melting pt: 77-78°C Solubility: Ethanol: 64 mg/mL warmed (199.97 mM); DMSO: 64 mg/mL (199.97 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kandul SV, et al. Fiziol Zh. 1986 Sep-Oct;32(5):579-84.

Note: This product is for in vitro research use only