

Phenothrin

Catalog No: #S1200


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phenothrin
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	350.45
Formulation	C23H26O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:CC(C)=CC1C(C(=O)OCC2=CC(OC3=CC=CC=C3)=CC=C2)C1(C)CPathways:OthersReceptor:OthersBoiling pt:437CMelting

pt:Solubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Nagy K, et al. Mutat Res Genet Toxicol Environ Mutagen. 2014 Aug;770:1-5.

Note: This product is for in vitro research use only