

Citicoline sodium

Catalog No: #S1208


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Citicoline sodium
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	510.31
Formulation	C ₁₄ H ₂₅ N ₄ NaO ₁₁ P ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:C[N+](C)(C)CCOP(=O)([O-])OP(=O)([O-])OC[C@@H]1[C@H]([C@H]([C@@H](O1)N2C=CC(=NC2=O)N)O)O.[Na+]Pathways:Others

Receptor:phosphatidylcholineBoiling pt:Melting pt:259-268C?Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Conant R, Schauss AG. Altern Med Rev. 2004 Mar;9(1):17-31.

Note: This product is for in vitro research use only