

Levosulpiride

Catalog No: #S1211


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Levosulpiride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Dopamine Receptor antagonist
Calculated MW	341.43
Formulation	C15H23N3O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: CCN1CCC[C@H]1CNC(=O)C2=C(C=CC(=C2)S(=O)(=O)N)OC
 Pathways: GPCR/G Protein;
 Neuroscience Receptor: Dopamine
 Boiling pt: Melting pt: 183-186°C
 Solubility: DMSO: 69 mg/mL (202.09 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Saporito F, et al. Int J Cardiol. 2014 Aug 1;175(2):383-4.

Note: This product is for in vitro research use only