

Ribostamycin sulfate

Catalog No: #S1214


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ribostamycin sulfate
Brief Description	Inhibitors
Purification	630µg/mg, BR
Target Name	Others
Calculated MW	552.55
Formulation	C17H34N4O10H2SO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C1[C@H]([C@@H]([C@H]([C@@H]([C@H]1N)O)[C@@H]2[C@@H]([C@H]([C@@H]([C@H](O2)CN)O)O)N)O[C@H]3[C@@H]([C@@H]([C@H]([C@H](O3)CO)O)O)O)N.OS(=O)(=O)O Pathways: Others Receptor: Others Boiling pt: 907.6C? Melting pt: 175-180 C Solubility: Water: 100 mg/mL (180.97 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Horibe T, et al. Biochem Biophys Res Commun, 2001, 289(5), 967-972.

Note: This product is for in vitro research use only