

Diethylstilbestrol

Catalog No: #S1216


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Diethylstilbestrol
Brief Description	Inhibitors
Purification	98.67%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	268.35
Formulation	C18H20O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CC/C(=C(/CC)\C1=CC=C(C=C1)O)/C2=CC=C(C=C2)OPathways:Endocrinology/HormonesReceptor:ERBoiling pt:407.1CMelting pt:170-172CSolubility:Ethanol: 54 mg/mL (201.22 mM); DMSO: 54 mg/mL (201.22 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Greschik H, et al. J Biol Chem. 2004 Aug 6;279(32):33639-46.

Note: This product is for in vitro research use only