

Sisomicin Sulfate

Catalog No: #S1219


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sisomicin Sulfate
Brief Description	Inhibitors
Purification	580lu/mg
Target Name	Others
Calculated MW	1385.43
Formulation	2C19H37N5O75H2SO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: OS(O)(=O)=O.CNC1C(O)C(OC2C(N)CC(N)C(OC3OC(CN)=CCC3N)C2O)OCC1(C)O Pathways: Others Receptor: protein-synthesis
 Boiling pt: 676.6°C Melting pt: 70-73°C Solubility: Water: 38mg/mL Appearance: Faint Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Hendriks JG, et al. Biomaterials. 2004 Feb;25(3):545-56.

Note: This product is for in vitro research use only