

## Cefuroxime sodium

Catalog No: #S1224


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cefuroxime sodium
Brief Description	Inhibitors
Purification	97.44%
Target Name	Others
Calculated MW	446.37
Formulation	C <sub>16</sub> H <sub>15</sub> N <sub>4</sub> NaO <sub>8</sub> S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: CO/N=C(/C1=CC=CO1)\C(=O)N[C@H]2[C@@H]3N(C2=O)C(=C(CS3)COC(=O)N)C(=O)[O-].[Na+] Pathways: Others Receptor: Others Boiling pt: Melting pt: 240-245°C Solubility: DMSO: 19mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Zhao L, et al. J Biomed Biotechnol. 2012;2012:507294.

Note: This product is for in vitro research use only