

## Pamabrom

Catalog No: #S1230


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pamabrom
Brief Description	Inhibitors
Purification	99.84%
Target Name	Others
Calculated MW	348.2
Formulation	C <sub>11</sub> H <sub>18</sub> BrN <sub>5</sub> O <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine system SMILES:CN1C2=C(NC(Br)=N2)C(=O)N(C)C1=O Pathways:Others Receptor:Others Boiling pt:Melting pt:>300C Solubility:Soluble in DMSO Appearance:White Powder Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Min JA, et al. Clin Exp Dermatol. 2009 Oct;34(7):e455-6.

Note: This product is for in vitro research use only