

Daidzein


Catalog No: #S1238

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Daidzein
Brief Description	Inhibitors
Purification	98.00%
Target Name	PPAR inhibitor
Calculated MW	254.24
Formulation	C15H10O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
SMILES: OC1=CC=C(C=C1)C1=COC2=C(C=CC(O)=C2)C1=O
Pathways: Metabolism
Receptor: PPAR
Boiling pt: Melting pt: 315-323°C
Solubility: DMSO: 51 mg/mL (200.59 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Feng G, et al. Int Immunopharmacol. 2015 Apr 15;26(2):392-400.

Note: This product is for in vitro research use only