

## Flumazenil

Catalog No: #S1240


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Flumazenil
Brief Description	Inhibitors
Purification	99.50%
Target Name	GABA Receptor antagonist
Calculated MW	303.29
Formulation	C15H14FN3O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous systemSMILES:CCOC(=O)C1=C2CN(C)C(=O)C3=C(C=CC(F)=C3)N2C=N1Pathways:Membrane transporter/Ion channel;  
NeuroscienceReceptor:GABARBoiling pt:528CMelting pt:201-203CSolubility:DMSO: 5 mg/mL (16.48 mM)Appearance:White PowderRemark:For  
obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C  
for several month.

## References

1. Wingrove PB, et al. Eur J Pharmacol. 2002 Feb 15;437(1-2):31-9.

Note: This product is for in vitro research use only