

Hydrocortisone acetate

Catalog No: #S1243


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hydrocortisone acetate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Annexin A; GR
Calculated MW	404.5
Formulation	C23H32O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:CC(=O)OCC(=O)[C@]1(CC[C@@H]2[C@@]1(C[C@@H]([C@H]3[C@H]2CCC4=CC(=O)CC[C@]34C)O)C)OP

Pathways:GPCR/G

Protein; MetabolismReceptor:Annexin A1; Glucocorticoid ReceptorBoiling pt:Melting pt:223CSolubility:Soluble in WaterAppearance:White

PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be

stored below -20 °C for several month.

References

1. Grossman R, et al. Ann N Y Acad Sci. 2006 Jul;1071:410-21.

Note: This product is for in vitro research use only