

Celiprolol hydrochloride

Catalog No: #S1244


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Celiprolol hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor antagonist
Calculated MW	415.96
Formulation	C ₂₀ H ₃₃ N ₃ O ₄ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System SMILES: Cl.CCN(CC)C(=O)NC1=CC(C(=O))=C(OCC(O)CNC(C)(C)C)C=C1 Pathways: GPCR/G

Protein Receptor: β -adrenergic receptor Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Abdelkrim MA, et al. Can J Physiol Pharmacol. 2013 Oct;91(10):791-6.

Note: This product is for in vitro research use only