

## Bornyl acetate

Catalog No: #S1246


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Bornyl acetate
Brief Description	Inhibitors
Purification	98.00%
Target Name	IL Receptor inhibitor
Calculated MW	196.29
Formulation	C <sub>12</sub> H <sub>20</sub> O <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CC(=O)OC1CC2CCC1(C)C2(C)C  
Pathways: Immunology/Inflammation  
Receptor: IL receptor  
Boiling pt: 228-231°C  
Melting pt: 27°C  
Solubility: DMSO: 10 mM  
Appearance: Liquid  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Yang H, et al. IUBMB Life. 2014 Dec;66(12):854-9.

Note: This product is for in vitro research use only