

Sulfanitran

Catalog No: #S1248


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sulfanitran
Brief Description	Inhibitors
Purification	98.30%
Target Name	Autophagy antagonist
Calculated MW	335.34
Formulation	C14H13N3O5S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: CC(=O)NC1=CC=C(C=C1)S(=O)(=O)NC1=CC=C(C=C1)[N+](=[O-])=O
 Pathways: Autophagy Receptor: DHPS
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Henry RJ, et al. Bacteriol Rev. 1943 Dec;7(4):175-262.

Note: This product is for in vitro research use only