

Beclamide

Catalog No: #S1249


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beclamide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	197.67
Formulation	C10H12ClNO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: C1CCC(=O)NCC1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: 386.9°C
Melting pt: Solubility: DMSO: 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Darmani NA, et al. J Pharm Pharmacol. 1988 Dec; 40(12):891-3.

Note: This product is for in vitro research use only