

## PAMBA


Catalog No: #S1252

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PAMBA
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	151.16
Formulation	C <sub>8</sub> H <sub>9</sub> NO <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular System  
SMILES: NCC1=CC=C(C(=O)O)C=C1  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: 300°C  
Solubility: Soluble in Water  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Ambasta SS, et al. Indian J Med Res. 1988 Jul;88:57-9.

Note: This product is for in vitro research use only