

Chlorindanol


Catalog No: #S1253

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chlorindanol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	168.62
Formulation	C ₉ H ₉ ClO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: OC1=C2CCCC2=C(Cl)C=C1
Pathways: Others
Receptor: Others
Boiling pt: 290°C
Melting pt: 91-93°C
Solubility: Soluble in DMSO
Appearance: Gray Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. BAILEY JH, et al. J Am Pharm Assoc Am Pharm Assoc (Baltim). 1959 Apr;48(4):212-6.

Note: This product is for in vitro research use only