

Ubenimex

Catalog No: #S1257


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ubenimex
Brief Description	Inhibitors
Purification	96.99%
Target Name	LTR inhibitor; Others inhibitor
Calculated MW	308.37
Formulation	C16H24N2O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CC(C)C[C@@H](C(=O)O)NC(=O)[C@H]([C@@H](CC1=CC=CC=C1)N)O
 Pathways: Immunology/Inflammation;
 Others Receptor: LTA(4)H-h; APN
 Boiling pt: 604.7°C
 Melting pt: 245°C
 Solubility: DMSO: 0.4 mg/mL (1.29 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sekine K, et al. Leukemia, 1999, 13(5), 729-734.

Note: This product is for in vitro research use only