

# Pyrantel pamoate


Catalog No: #S1259

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pyrantel pamoate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	594.69
Formulation	C34H30N2O6S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

### Research

Area: Immunology SMILES: CN1CCCN=C1/C=C/C2=CC=CS2.C1=CC=C2C(=C1)C=C(C(=C2CC3=C(C(=CC4=CC=CC=C43)C(=O)O)O)O)C(=O)OPat  
 hways: Others Receptor: Antiparasitic Boiling pt: 266-267°C Melting pt: Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Salman AB. J Pediatr Surg. 1997 Apr;32(4): 585-7.

Note: This product is for in vitro research use only