

Farlutin

Catalog No: #S1261


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Farlutin
Brief Description	Inhibitors
Purification	99.45%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	386.54
Formulation	C24H34O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@H]1C[C@@H]2[C@H](CC[C@]3([C@H]2CC[C@@]3(C(=O)C)OC(=O)C)[C@@]4(C1=CC(=O)CC4)CPathways:Endocrinology/HormonesReceptor:ER; PRBoiling pt:496.4CMelting pt:206-207CSolubility:Ethanol: 12 mg/mL (31.04 mM); DMSO: 12 mg/mL (31.04 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Risch HA, et al. Y Epidemiol Biomarkers Prev. 2006 Sep;15(9):1738-41.

Note: This product is for in vitro research use only