

Ciprofibrate

Catalog No: #S1264


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ciprofibrate
Brief Description	Inhibitors
Purification	99.50%
Target Name	PPAR agonist
Calculated MW	289.16
Formulation	C ₁₃ H ₁₄ Cl ₂ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular system
 SMILES: CC(C)(C(=O)O)OC1=CC=C(C=C1)C2CC2(Cl)Cl
 Pathways: Metabolism
 Receptor: PPAR
 Boiling pt: 424.9°C
 Melting pt: 114-118°C
 Solubility: Soluble in DMSO
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Domingos H, et al. Arq Bras Cardiol. 2012 Nov;99(5):997-1007.

Note: This product is for in vitro research use only