

Dixanthogen


Catalog No: #S1268

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dixanthogen
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	242.4
Formulation	C6H10O2S4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: CCOC(=S)SSC(=S)OCC Pathways: Others Receptor: Others Boiling pt: 134°C Melting pt: Solubility: Soluble in DMSO Appearance: Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. DANBOLT N. Acta Derm Venereol. 1945; 26(1):103.

Note: This product is for in vitro research use only