

Diacerein

Catalog No: #S1273


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Diacerein
Brief Description	Inhibitors
Purification	98.00%
Target Name	IL Receptor inhibitor
Calculated MW	368.29
Formulation	C19H12O8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous

system SMILES: CC(=O)OC1=CC=CC2=C1C(=O)C1=C(OC(C)=O)C=C(C=C1C2=O)C(O)=O Pathways: Immunology/Inflammation Receptor: IL-1 β Binding pt: Melting pt: 217-218°C Solubility: DMSO: 48 mg/mL (130.33 mM) Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Yaron M, et al. Osteoarthritis Cartilage, 1999, 7(3), 272-280.

Note: This product is for in vitro research use only