

## Benzbromarone

Catalog No: #S1274


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Benzbromarone
Brief Description	Inhibitors
Purification	98.00%
Target Name	P450 inhibitor
Calculated MW	424.08
Formulation	C17H12Br2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism

system SMILES: CCC1=C(C(=O)C2=CC(Br)=C(O)C(Br)=C2)C2=CC=CC=C2O1 Pathways: Metabolism Receptor: CYP2C9

Boling pt: Melting pt: Solubility: Ethanol: 9 mg/mL (21.22 mM); DMSO: 85 mg/mL (200.43 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. McDonald MG, et al. Chem Res Toxicol, 2007, 20(12), 1833-18342.

Note: This product is for in vitro research use only