

Carteolol hydrochloride

Catalog No: #S1282


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Carteolol hydrochloride
Brief Description	Inhibitors
Purification	99.91%
Target Name	Adrenergic Receptor antagonist
Calculated MW	328.83
Formulation	C ₁₆ H ₂₄ N ₂ O ₃ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System SMILES: Cl.CC(C)(C)NCC(O)COC1=CC=CC2=C1CCC(=O)N2 Pathways: GPCR/G

Protein Receptor: β -adrenergic receptor Boiling pt: 518.6°C Melting pt: 278°C Solubility: Water: 65 mg/mL (197.67 mM); DMSO: 65 mg/mL (197.67

mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Wheeldon NM, et al. Br J Clin Pharmacol. 1992 Apr;33(4):411-6.

Note: This product is for in vitro research use only