

## Anethole trithione

Catalog No: #S1285


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Anethole trithione
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	240.35
Formulation	C10H8OS3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system  
 SMILES: CC=CC1=CC=C(OC=S)C(=S)C1=S  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: 398.1°C  
 Melting pt: [not specified]  
 Solubility: DMSO: 10 mM  
 Appearance: Red Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Ukai Y, et al. *Arzneimittelforschung*. 1988 Oct; 38(10):1460-5.

Note: This product is for in vitro research use only