

Vincamine


Catalog No: #S1286

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vincamine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	354.44
Formulation	C ₂₁ H ₂₆ N ₂ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: CC[C@@]12CCCN3[C@@H]1C4=C(CC3)C5=CC=CC=C5N4[C@](C2)(C(=O)OC)O Pathways: Others Receptor: Others Boiling pt: 508.9°C Melting pt: 232°C Solubility: Water: 62 mg/L Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Pu H, et al. Luminescence. 2014 Aug;29(5):471-9.

Note: This product is for in vitro research use only