

TTNPB (Arotinoid Acid)


Catalog No: #S1288

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TTNPB (Arotinoid Acid)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Retinoid Receptor agonist
Calculated MW	348.48
Formulation	C ₂₄ H ₂₈ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C/C(=C/C1=CC=C(C(=C1)C(=O)O)/C2=CC3=C(C=C2)C(CCC3(C)C)(C)C
Pathways: Metabolism
Receptor: RAR α ; RAR β ; RAR γ ; RAR γ
Boiling pt: Melting pt: Solubility: DMSO: 15 mg/mL warmed (43.04 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Pignatello MA, et al. Toxicol Appl Pharmacol. 1997, 142(2), 319-327.

Note: This product is for in vitro research use only