

Eugenol


Catalog No: #S1297

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Eugenol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	164.2
Formulation	C10H12O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC1=C(O)C=CC(CC=C)=C1
Pathways: Others
Receptor: Others
Boiling pt: 254°C
Melting pt: 12-10°C
Solubility: DMSO: 10 mM
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Jadhav BK, et al. Drug Dev Ind Pharm. 2004 Feb;30(2):195-203.

Note: This product is for in vitro research use only