

Clofibrate


Catalog No: #S1298

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Clofibrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	PPAR agonist
Calculated MW	242.7
Formulation	C ₁₂ H ₁₅ ClO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System SMILES: CCOC(=O)C(C)(C)OC1=CC=C(Cl)C=C1 Pathways: Metabolism Receptor: PPAR α Boiling pt: 274.8°C Melting pt: Solubility: DMSO: 7 mg/mL; Water:

References

1. Casas F, et al. FEBS Lett. 2000, 478(1-2):4-8.

Note: This product is for in vitro research use only