

## Roniacol

Catalog No: #S1299


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Roniacol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	109.13
Formulation	C6H7NO
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism systemSMILES:OCC1=CN=CC=C1Pathways:OthersReceptor:OthersBoiling pt:154CMelting pt:-7CSolubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Dujovne CA, et al. Am J Med Sci. 1979 May-Jun;277(3):255-61.

Note: This product is for in vitro research use only