

Propofol

Catalog No: #S1300


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Propofol
Brief Description	Inhibitors
Purification	99.13%
Target Name	GABA Receptor; Sodium Channel inhibitor
Calculated MW	178.28
Formulation	C12H18O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system SMILES: CC(C)C1=CC=CC(C(C)C)=C1O Pathways: Membrane transporter/Ion channel;

Neuroscience Receptor: GABAR; Sodium Channel Boiling pt: 256°C Melting pt: 18°C Solubility: Soluble in Water Appearance: Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Franks NP. Br J Pharmacol. 2006 Jan;147 Suppl 1:S72-81.

Note: This product is for in vitro research use only