

Guaiacol


Catalog No: #S1301

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Guaiacol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	124.14
Formulation	C7H8O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory systemSMILES:COC1=CC=CC=C1OPathways:OthersReceptor:OthersBoiling pt:205CMelting pt:26-29CSolubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Dorfner R, et al. J Agric Food Chem. 2003 Sep 10;51(19):5768-73.

Note: This product is for in vitro research use only