

## Auranofin

Catalog No: #S1303


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Auranofin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	678.48
Formulation	C <sub>20</sub> H <sub>34</sub> AuO <sub>9</sub> PS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: [Au+].CCP(CC)CC.CC(=O)OCC1OC([S-])C(OC(C)=O)C(OC(C)=O)C1OC(C)=O Pathways: Others Receptor: TrxRB  
 Boiling pt: 425.5°C Melting pt: 103-105°C Solubility: DMSO: 10 mM Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Rigobello MP, et al. Free Radic Res. 2005 Jul; 39(7):687-95.

Note: This product is for in vitro research use only