

Deracil

Catalog No: #S1310


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Deracil
Brief Description	Inhibitors
Purification	98.00%
Target Name	Thyroid hormone receptor(THR) inhibitor
Calculated MW	128.16
Formulation	C ₄ H ₄ N ₂ O ₂ S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:O=C1NC(=S)NC=C1Pathways:Endocrinology/HormonesReceptor:Thyroid peripheral conversionBoiling pt:337.2CMelting pt:>300CSolubility:DMSO: 26 mg/mL (202.88 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Palumbo A, et al. FEBS Lett, 2000, 485(2-3), 109-112.

Note: This product is for in vitro research use only