

4-aminobenzoic acid


Catalog No: #S1311

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-aminobenzoic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	137.14
Formulation	C7H7NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: NC1=CC=C(C=C1)C(=O)O
Pathways: Others
Receptor: p-HBHB
Boiling pt: Melting pt: 187-189°C
Solubility: soluble in Water
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Benali-Cherif R, et al. Acta Crystallogr C Struct Chem. 2014 Mar;70(Pt 3):323-5.

Note: This product is for in vitro research use only