

## Cinchonidine

Catalog No: #S1317


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cinchonidine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	294.39
Formulation	C19H22N2O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: C=C[C@H]1CN2CC[C@H]1C[C@H]2[C@@H](C3=CC=NC4=CC=CC=C34)O Pathways: Others Receptor: Others Boiling pt: 464.5°C Melting pt: 204-206°C Solubility: Ethanol: 59 mg/mL (200.41 mM); DMSO: 59 mg/mL (200.41 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Fukiya K, et al. Drug Metab Dispos. 2010 Feb; 38(2):302-7.

Note: This product is for in vitro research use only