

Diatrizoate sodium

Catalog No: #S1319


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Diatrizoate sodium
Brief Description	Inhibitors
Purification	99.89%
Target Name	Others
Calculated MW	635.9
Formulation	C ₁₁ H ₈ I ₃ N ₂ NaO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: [Na+].CC(=O)NC1=C(I)C(C([O-])=O)=C(I)C(NC(C)=O)=C1
 Pathways: Others
 Receptor: Others
 Boiling pt: 614.1
 Melting pt: Solubility: Water: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Abd El-Rahman MK, et al. Spectrochim Acta A Mol Biomol Spectrosc. 2015 Feb 5;136 Pt B:1167-74.

Note: This product is for in vitro research use only