

Bendazol

Catalog No: #S1322


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bendazol
Brief Description	Inhibitors
Purification	98.70%
Target Name	Dopamine Receptor antagonist
Calculated MW	208.26
Formulation	C ₁₄ H ₁₂ N ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system SMILES: C(C1=NC2=CC=CC=C2N1)C1=CC=CC=C1 Pathways: GPCR/G Protein;

Neuroscience Receptor: Dopamine Boiling pt: 453.5°C Melting pt: 187°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kuzubova EA, et al. Eksp Klin Farmakol. 2007 Mar-Apr;70(2):37-9.

Note: This product is for in vitro research use only