

## Flavin mononucleotide

Catalog No: #S1331


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Flavin mononucleotide
Brief Description	Inhibitors
Purification	92.15%
Target Name	Vitamin inhibitor
Calculated MW	478.33
Formulation	C17H20N4NaO9P
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: CC1=CC2=C(C=C1C)N(C3=NC(=O)NC(=O)C3=N2)C[C@@H]([C@@H]([C@@H](COP(=O)(O)O)O)O)O.[Na+] Pathways: Metabolism Receptor: VB2 Boiling pt: Melting pt: >300°C Solubility: Water: >5mM Appearance: Orange Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Shen K, Huang XE. Asian Pac J Y Prev. 2015;16(4):1525-7.

Note: This product is for in vitro research use only