

Pamidronate disodium salt


Catalog No: #S1337

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pamidronate disodium salt
Brief Description	Inhibitors
Purification	99.93%
Target Name	HBV; Others
Calculated MW	279.03
Formulation	C3H9NNa2O7P2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: O.O.O.O.O.O.O.[Na+].[Na+].NCCC(O)(P(O)([O-])=O)P(O)([O-])=O
 Pathways: Microbiology&virology;
 Others Receptor: HAP; FPPS
 Boiling pt: Melting pt: 300C
 Solubility: Water: 56 mg/mL (200.69 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Bergstrom JD, et al. Arch Biochem Biophys. 2000 Jan 1;373(1):231-41.

Note: This product is for in vitro research use only