

N-Carboethoxy-4-piperidone


Catalog No: #S1340

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-Carboethoxy-4-piperidone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	171.19
Formulation	C8H13NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCOC(=O)N1CCC(=O)CC1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Eric Gundersen et al. *Bioorganic & Medicinal Chemistry Letters*, 15(7), 1891-1894 (2005)

Note: This product is for in vitro research use only