

## Nortropine


Catalog No: #S1341

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nortropine
Brief Description	Inhibitors
Purification	99.00%
Target Name	AChR antagonist
Calculated MW	127.18
Formulation	C7H13NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: C1C[C@H]2CC(C[C@@H]1N2)O  
Pathways: Neuroscience  
Receptor: mAChR  
Boiling pt: Melting pt: Solubility: Soluble in DMSO  
Appearance: White Crystal  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. BA Bartholomew et al. Appl. Environ. Microbiol. September 1996

Note: This product is for in vitro research use only