

3-Chloropropiophenone


Catalog No: #S1342

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | 3-Chloropropiophenone |
| Brief Description | Inhibitors |
| Purification | 98%+ |
| Target Name | Others |
| Calculated MW | 168.62 |
| Formulation | C ₉ H ₉ ClO |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Other
SMILES: C1CCC(=O)C1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Gen-Sheng Y, et al. *Biotechnol Lett.* 2009 Dec;31(12):1879-83.

Note: This product is for in vitro research use only