

Leucocrystal Violet


Catalog No: #S1343

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Leucocrystal Violet
Brief Description	Inhibitors
Purification	99%+
Target Name	Others
Calculated MW	373.53
Formulation	C25H31N3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Other
 SMILES:CN(C)C1=CC=C(C=C1)C(C1=CC=C(C=C1)N(C)C)C1=CC=C(C=C1)N(C)C
 Pathways:Others
 Receptor:Others
 Boiling pt:
 Melting pt:
 Solubility:Soluble in Water
 Appearance:White Crystal
 Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several month.

References

1. Ju S, et al. Food Chem. 2015 Oct 15;185:479-87.

Note: This product is for in vitro research use only