

4-Chloro-4'-hydroxybenzophenone


Catalog No: #S1344

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Chloro-4'-hydroxybenzophenone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	232.66
Formulation	C13H9ClO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC1=CC=C(C(=O)C1=CC=C(Cl)C=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Cox PJ, et al. Acta Crystallogr B. 2008 Apr;64(Pt 2):206-16.

Note: This product is for in vitro research use only