

1-Phenyl-1,2-ethanediol

Catalog No: #S1346


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	1-Phenyl-1,2-ethanediol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	138.16
Formulation	C8H10O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:OCC(O)C1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited
solubilityAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a
while.Stock solution can be stored below -20 °C for several month.

References

- 1.Diana Poli et al. Chemical Research in Toxicology, 17(1), 104-109 (2004)

Note: This product is for in vitro research use only