

4-Butoxyphenol

Catalog No: #S1350


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	4-Butoxyphenol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	166.22
Formulation	C10H14O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CCCCOC1=CC=C(O)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited
solubilityAppearance:Brown Sheet CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath
for a while.Stock solution can be stored below -20 °C for several month.

References

1. Davydov A, Grönlund A. Biochem Biophys Res Commun. 1999 May 10;258(2):322-5.

Note: This product is for in vitro research use only