

DL-Mevalonolactone


Catalog No: #S1353

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | DL-Mevalonolactone |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Others |
| Calculated MW | 130.14 |
| Formulation | C6H10O3 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Other
SMILES: CC1(O)CCOC(=O)C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Endo A, et al. J Lipid Res. 1992 Nov; 33(11):1569-82.

Note: This product is for in vitro research use only